

Data Sheet

Mixing / Magnetic Stirrers



topolino

Extremely convenient magnetic mini-stirrer for mixing quantities up to 250 ml.

- Durable, brushless motor
- Infinitely variable speed range
- High magnetic adhesion

Accessories: TRIKA® 25 Magnetic stirring bar, RSE Stirring bar remover, IKAFLON® 10 Magnetic stirring bar, IKAFLON® 30 Magnetic stirring bar, IKAFLON® 20 Magnetic stirring bar, IKAFLON® 25 Magnetic stirring bar, IKAFLON® 15 Magnetic stirring bar

Technical Data

| | |
|---|----------------|
| Number of stirring positions | 1 |
| Stirring quantity max. per stirring position (H ₂ O) [l] | 0.25 |
| Stirring quantity max. (H ₂ O) [l] | 0.25 |
| Motor rating input [W] | 1 |
| Motor rating output [W] | 0.8 |
| Speed range [rpm] | 300 - 1800 |
| Stirring bar length max. [mm] | 30 |
| Speed control | stepless |
| Set-up plate material | synthetic (PP) |
| Set-up plate dimensions [mm] | Ø 80 |
| Dimensions (W x H x D) [mm] | 95 x 37 x 115 |
| Weight [kg] | 0.32 |
| Permissible ambient temperature [°C] | 5 - 40 |
| Permissible relative moisture [%] | 80 |
| Protection class according to DIN EN 60529 | IP 21 |
| RS 232 interface | no |
| Analog output | no |
| Voltage [V] | 100 - 240 |
| Frequency [Hz] | 50/60 |
| Power input [W] | 10 |
| Ident. No. | 3368000 |