

electronic pipette controller

Profiller™ *electro* 446


Pipette controller fitting any hand in a very comfortable way. Shape is ideal for working with glass or plastic pipettes within laminar flow or microbiological safety cabinet. A powerful, yet quiet motor and a switch to set both aspiration and dispensing speeds enable maximal working efficiency even with large pipetting volumes. One-year warranty.

Product features

- Powerful pump - fills 25 mL in 3 seconds
- Simple, fast operation
- Comfortable finger triggers
- Low battery light
- NiMH battery rechargeable during use
- Interchangeable protection membrane filter
- Stand suitable for bench or wall mounting



Pipette volume up to 100 mL



Finger triggers



Variable speed



Battery charge

Convenient stand

Suitable for bench or wall mounting, stand houses instrument with or without a pipette in place.



Ordering information - Instrument

Supplied with charger, two spare hydrophobic filters (0.45µm and 0.2µm), bench/wall stand and instruction manual.

Description	Pack.	Cat. No.
Profiller™ <i>electro</i> 446	1 / pk	446.100.x

Replace x by country code for correct plug type depending on user location : E = Europe, G = UK, U = USA-Japan

Ordering information - Accessories

Description	Pack.	Cat. No.
Spare bench stand for Profiller™ 446	1 / pk	320.446
Replacement 1.2V NiMH battery	1 / pk	900.916
Fast charger + cord, Europe style plug	1 / pk	900.946E
Fast charger + cord, UK style plug	1 / pk	900.946G
Fast charger + cord, USA-Japan style	1 / pk	900.946U
PTFE Membrane filter, 0.45 µm	5 / pk	322.435
PTFE Membrane filter, 0.2 µm	5 / pk	322.400
Spare silicone pipette holder	1 / pk	1.435.01
Spare nose cone white	3 / pk	1.435.02
Alt. nose cones pink, violet, green (1 ea.)	3 / pk	1.435.03

For more information or for an overview of the complete Socorex product range, visit

www.socorex.com