

manual pipette controllers

Profiller™ manual 435 / 436

Robust and easy to use, both models offer ergonomic shape and light weight. They provide smooth control in both aspirating and dispensing, using glass or plastic pipettes of volume up to 100 mL. One-year warranty.

Product features

- Squeezable large capacity bulb
- Soft action lever button for precise control
- Efficient blow-out system
- PTFE protection membrane filter
- Maintenance free
- Alternative colour coded nose cones



Pipette volume up to 100 mL



Squeezable large capacity bulbs



Precise control



Liquid blow-out



Protection filter



Colour cones

Ordering information - Instruments

Description	Packaging	Cat. No.
Profiller™ manual 435, egg style grey bulb	1 / pk	435.100
Profiller™ manual 436, flat shaped blue bulb	1 / pk	436.100

Ordering information - Accessories

Description	Packaging	Cat. No.
PTFE Membrane filter, 0.45 µm	5 / pk	322.435
PTFE Membrane filter, 0.2 µm	5 / pk	322.400
Spare silicone pipette holder	1 / pk	1.435.01
Spare nose cone in white	3 / pk	1.435.02
Alternative nose cones in pink, violet and green (1 each)	3 / pk	1.435.03