

BABY SCALE

PRODUCT NUMBER 8320

Medicine, health and wellness



- › Foldable to half its length
- › Integrated carrier handle
- › Lightweight and easy to handle
- › Eight-point measuring system guarantees exact weight evaluation

BABY SCALE

PRODUCT NUMBER 8320



Features

The first baby scale with practical collapsible function for parents, midwives and pharmacists can easily reduce its length by half and is impressively light. Furthermore, the integrated carrier handle and successful beauty case design, make it especially manageable and easy to transport. Even with an unsettled baby, the eight-point measuring system guarantees exact weight evaluation, easily legible on a large numeric display. The baby scale contains an environmentally-friendly auto-off function, leading to longer battery life.

Technical Data

- › Electronic weighing system
- › Eight load cells
- › Switchover to lb/oz units
- › Power supply: 2 x 1.5V AA cells
- › Scale switches off automatically for longer battery life

Design

- › Tare function
- › Steady weight display through filtering and hold function
- › Large 20mm LCD digits for easy reading
- › Small space requirement and easy to transport when folded
- › Folded scale held together by magnets
- › Easy and fast opening: Open with carrier handles and remove the base
- › Plastic housing
- › Patent pending

Ordering information

Order Number	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
8320.01.001	15	10	375 x 315 x 105 (closed) 740 x 315 x 53 (open)	1,9

Your Soehnle Professional dealer:

www.soehnle-professional.com

BABY SCALE

PRODUCT NUMBER 8310

Medicine, health and wellness



- > Weighing range 20kg
- > Division: 10g
- > Hold function
- > Automatic power-off to save battery capacity

BABY SCALE

PRODUCT NUMBER 8310

Features

- › Beautifully designed handy baby scale made of easy to clean plastic
- › Removable baby weighing tray: can be used as personal scale for weighing children up to 20kg.
- › Quattro-Tronic: fully electronic scale with 4 load cells for rapid precise weight measurement
- › No wear parts and maintenance-free
- › Easy to read LCD with 20mm high digits
- › Simple one-key operation for on/off and tare
- › Weight lock function: weight display remains stable if baby moves about.
- › Blankets can be tared.
- › Lightweight sturdy carry case to protect the scale when transported
- › Battery power supply with 9V block battery and energy saving automatic power-off

Ordering information

Order Numbers	Order Numbers	Weighing Range	Divisions	Dimensions	Intrinsic weight
	Carry case	kg	g	mm	kg
8310.00.090	5030.00.071	20	10	530 X 240 X 100	4

Your Soehnle Professional dealer:

www.soehnle-professional.com

BABY SCALE

PRODUCT NUMBER 7726.03

Medicine, health and wellness



- > Approval Class III
- > MPG Class 1
- > Backlit LCD with 20mm high digits
- > Hold function

BABY SCALE

PRODUCT NUMBER 7726.03



Features

Two-range scale with high precision and reliability, stable weight display by damping and hold function, high stability baby weighing tray made of easy to clean plastic.

- › Easy to operate
- › High contrast, backlit LCD with 20mm high digits
- › Automatic zeroing on power-on
- › Tare button (to tare blankets, etc.)
- › Print and data transfer button
- › RS 232 data interface (V24)
- › Mains power supply or rechargeable battery
- › With weight lock function
- › Unit indicator switchover kg/lb

Accessories

- › Baby height rod (Order No. 5001.10.001)

Ordering information

Order Numbers	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7726.03.001	15	up to 6kg = 2g up 6kg = 5g	530 x 320 x 180	6.6

Your Soehnle Professional dealer:

www.soehnle-professional.com