

# ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7702

Medicine, health and wellness



- › Electronic chair weigher
- › Approved to Class III, MPG Class 1
- › Hinged footrests and armrests
- › No wear parts and maintenance-free
- › Weighing range 200kg
- › BMI and dialysis functions

# ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7702



## Features

Electronic chair weighers to meet very high requirements – besides approval to EU Class III and to MPG, the main features include high precision and reliability, user-friendliness and easy readability. These features are ideally suited for hospitals, old people's homes and nursing homes for weighing patients who have difficulty walking or standing.

### Design

- > Approved to EU Class III, MPG Class 1
- > Seat made of easy to clean plastic, ergonomically shaped
- > Hinged footrests and armrests
- > 4 wheels, of which 2 lockable steering wheels
- > Digital indicator: housing made of aluminium sections with plastic parts
- > 7-digit LCD, backlit, digit height 20mm
- > Two-range version
- > Units in kg or lb

### Technical data

- > Mains power: power adapter input 230VAC, 50Hz, output 9VDC, 300mA
- > Battery power supply for approx. 160 hrs by 4 x 1.5V Mignon cells
- > Auto-Off function (time adjustable)
- > Dimensions: lxxh 580 x 1100 x 930mm
- > Unit switchover kg/lb

### Options

- > Charge adapter with rechargeable batteries for approx. 120 hrs service time (Order No. 2557.04.001)
- > Pouch for power adapter for attachment to wheelchair scale (Order No. 2507.03.001)
- > RS 232 serial interface for connection to PC/printer

## Ordering information

Order Numbers	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7702.01.001 (without interface)	200	bis 150kg = 50g ab 150kg = 100g	580 x 1100 x 930	27
7702.01.002 (with interface)	200	bis 150kg = 50g ab 150kg = 100g	580 x 1100 x 930	27

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)

# ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7802

Medicine, health and wellness



- › Weighing range 200kg
- › Hinged footrests and armrests
- › No wear parts and maintenance-free

# ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7802



## Features

The mobile Electronic Chair Weigher 7802 is the ideal scale for checking the weight of people who have difficulty walking or standing in areas not subject to the Medical Product Law. The Soehnle Professional electronic chair weigher is easy to use and has a sturdy design, large indicator, comfortable seat and ergonomic timeless design.

### Design

- › Seat made of easy to clean plastic, ergonomically shaped
- › Hinged footrests
- › Hinged armrests
- › 4 wheels, of which 2 lockable steering wheels
- › Digital indicator: housing made of impact resistant plastic
- › No wear parts and maintenance-free
- › Simple one-key operation

### Technical data

- › Mains power data: 230V, 50/60Hz, 10mA, 2.2VA
- › Battery power supply for approx. 150 hrs with 4 x 1.5V Mignon cells
- › Auto-Off function in battery power mode
- › Dimensions: lxxh 580 x 1100 x 930mm
- › Unit switchover kg/lb

### Options

- › Charge adapter with rechargeable batteries for approx. 80 hrs continuous service time (Order No. 2557.03.001)

## Ordering information

Order Numbers	Order Numbers	Weighing Range	Division	Dimensions	Weight
rechargeable battery	rechargeable battery	kg	g	mm	kg
mains power supply	incl. charge adapter				
7802.01.001	7802.01.002	200	100	580 x 1100 x 930	27
incl. power supply	incl. 2557.03.001				

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)